



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Yeast Autolysate Supplement

Product Code: MS2027

Application: An enrichment supplement recommended for the cultivation of *Neisseria* species.

Composition**

Per vial sufficient for 500 ml medium

Ingredients	Concentration
Yeast autolysate	5g
Dextrose	0.500g
Sodium bicarbonate	0.075g

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well & aseptically add to sterile, molten, cooled (45-50°C) 490 ml of GC Agar Base DM1434 / GC MiVeg Agar Base VM1434 containing Haemoglobin MS2022 along with GC supplement w/ antibiotics MS2021 or In Chocolate Agar Base DM1103 / Chocolate MiVeg Agar Base VM1103 along with 5% v/v defibrinated blood (after addition of blood mix and heat at 80°C for 10mins) or In 320 ml of NYC Agar Base DM2348 alongwith 100ml sedimented horse blood cells, 60ml citrated Horse plasma and NYC supplement MS2150 . Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.